

**American College of Physicians - Internal Medicine Meeting 2026
San Francisco, CA**

POCUS for LV Function

Faculty Information

Director:

Woo J. Moon, DO, FACP

Shane Robert Kinard, MD, FACP

Damien Ricklis, MD, FACP

Ian R. Ross, MD, FACP

Hillary Kieran Deveaux, MBBS, Resident/Fellow Member

Posted Date: March 10, 2026

©2026 American College of Physicians. All rights reserved. Reproduction of presentations, or print or electronic material associated with presentations, is prohibited without written permission from the ACP.

Any use of program content, the name of a speaker and/or program title, or the name of ACP without the written consent of ACP is prohibited. For purposes of the preceding sentence, "program content" includes, but is not limited to, oral presentations, audiovisual materials used by speakers, program handouts, and/or summaries of the same. This rule applies before, after, and during the activity

POCUS for Left Ventricular Systolic Function

Shane R Kinard, MD, FACP. *Penn State College of Medicine*

Woo J Moon, DO, FACP. *Saint Louis University School of Medicine*

Damien L Ricklis, MD, FACP. *University of Kansas School of Medicine*

Ian R Ross, MD, FACP. *Washington University School of Medicine in St. Louis*



1

Disclosures

Visit any speaker's profile within the *ACP Meeting* mobile app or the meeting's web platform to view disclosure of relevant financial relationships.



2

Learning Objectives

- Describe standard cardiac ultrasound image orientation
- Acquire parasternal long axis, parasternal short axis, apical four chamber and subcostal four chamber views
- Estimate left ventricular systolic function using visual estimation and E-point septal separation (EPSS)
- Recognize reduced left ventricular systolic function on POCUS imaging

3


Probe Selection

- Phased-array probe (1-5 MHz)
 - Low frequency for adequate depth
- Small footprint for intercostal imaging
- High frame rate to evaluate cardiac motion



4

Probe Manipulation



The diagram illustrates four probe manipulation techniques with corresponding ultrasound images:

- Slide:** Shows two probe positions with a red arrow indicating lateral movement. The ultrasound image shows a vessel lumen with the probe marker on the left.
- Rotate:** Shows a probe being rotated around its axis, indicated by a fan of red arrows. The ultrasound image shows the vessel lumen with the probe marker on the left.
- Fan/Tilt:** Shows a probe being tilted to fan the beam. The ultrasound image shows a wider field of view of the vessel lumen with the probe marker on the left.
- Rock:** Shows a probe being rocked to change the view. The ultrasound image shows a different view of the vessel lumen with the probe marker on the left.

ACP American College of Physicians®
Leading Internal Medicine, Improving Lives

5

Cardiac Image Orientation

- Longitudinal views (i.e. PLAX): Probe marker cephalad
- Transverse views (i.e. PSAX): Probe marker leftward
- Screen indicator: Upper right corner
- Superficial structures appear at the top of the image

6

Cardiac Image Orientation


7

Cardiac POCUS Exam Views

- Parasternal long axis
- Parasternal short axis
- Apical four chamber
- Subcostal four chamber
- Inferior vena cava

8


Parasternal Long Axis Image Acquisition




- Probe position: Left parasternal border
- Marker: Toward the right shoulder
- Slide caudad until the heart is visualized
- Rock and fan to optimize the image
- Adjust depth to include the descending thoracic aorta

9

Parasternal Long Axis Image Acquisition




P5-1
Focused Cardiac
MI: 1.1 TIS: 0.8



2D: G: 28
Pen DR: 0
THI

10



Parasternal Short Axis Image Acquisition



- Probe: Rotate 90° clockwise from parasternal long axis
- Marker: Toward the left shoulder
- Fan or slide to obtain the mid-ventricular papillary muscle level

11

Parasternal Short Axis Image Acquisition



P5-1
Cardiac
MI: 1.1 TIS: 0.8

2D: G: 50
Pen DR: 0
THI

12

Left Ventricular Systolic Function Assessment

Left Ventricular Systolic Function


- Assessed qualitatively: Hyperdynamic → Severely reduced
- Based on three criteria:
 - 1) Endocardial excursion
 - 2) Myocardial thickening
 - 3) Anterior mitral valve leaflet motion

Left Ventricular Systolic Function

1. Endocardial excursion:

Symmetric endocardial motion towards the LV chamber

- Parasternal Long Axis
- Parasternal Short Axis



P5-1
Cardiac
MI: 0.9 TIS: 0.8


2D: G: 10
Pen DR: 0
THI

Left Ventricular Systolic Function

1. Endocardial excursion:

Symmetric endocardial motion towards the LV chamber

- Parasternal Long Axis
- Parasternal Short Axis



P5-1
Focused Cardiac
MI: 0.9 TIS: 0.8


2D: G: 42
Pen DR: 0
THI

Left Ventricular Systolic Function

2. Myocardial thickening:

>40% increase in thickness during systole

- Parasternal Long Axis
- Parasternal Short Axis



P5-1
Cardiac
MI: 0.9 TIS: 0.8


2D: G: 10
Pen DR: 0
THI

Left Ventricular Systolic Function

2. Myocardial thickening:

>40% increase in thickness during systole

- Parasternal Long Axis
- Parasternal Short Axis



P5-1
Focused Cardiac
MI: 0.9 TIS: 0.8


2D: G: 42
Pen DR: 0
THI

Left Ventricular Systolic Function

3. Anterior mitral valve leaflet motion:

< 1cm separation between the mitral valve tip and septum (E-point septal separation)

- Parasternal Long Axis




P5-1
Focused Cardiac
MI: 1.1 TIS: 0.8

2D: G: 28
Pen DR: 0
THI

M-Mode


- **What:** Evaluates motion of cardiac structures over time
- **How:** Obtain a B-mode image, then place the M-mode cursor along the structure of interest
- **Why:** High temporal resolution, ideal for rapidly moving structures
- Y-axis = depth, X-axis = time, visually depicts motion as depth vs. time

E-Point Septal Separation: M-Mode



21

E-Point Septal Separation: M-Mode



22

Take Home Points

- Use standard cardiac convention with indicator in the upper right of the screen
- Visually estimate function based on endocardial excursion and myocardial wall thickening
- E-point septal separation ≥ 10 mm suggests reduced systolic function
 - +LR 9.8, -LR 0.2 (Núñez-Ramos JA et al, POCUS J. 2022;7(1): 160-165)
- Learn more: [ACP POCUS Module #3](#)


References

1. Kimura BJ, Shaw DJ, Amundson SA, Phan JN, Blanchard DG, DeMaria AN. Cardiac limited ultrasound examination techniques to augment the bedside cardiac physical examination. *J Ultrasound Med.* 2015;34(9):1683-1690. doi:10.7863/ultra.15.14.09002
2. Marbach JA, Almufleh A, Di Santo P, et al. Comparative Accuracy of Focused Cardiac Ultrasonography and Clinical Examination for Left Ventricular Dysfunction and Valvular Heart Disease: A Systematic Review and Meta-analysis. *Ann Intern Med.* 2019;171(4):264-272. doi:10.7326/M19-1337
3. D'Cruz IA, Lalimalani GG, Sambasivan V, Cohen HC, Glick G. The superiority of mitral E point-ventricular septum separation to other echocardiographic indicators of left ventricular performance. *Clin Cardiol.* 1979;2(2):140-145. doi:10.1002/clc.4960020210
4. Rownak Jahan Tamanna, Shabnam Jahan Hoque, Faisal Mohammed Pasha. E-Point Septal Separation: A Bedside Tool for Emergency Physician Assessment of Left Ventricular Ejection Fraction. *Bangladesh Crit Care J.* 2023;11(2):90-94. doi:10.3329/bccj.v11i2.69181
5. Soni MD MS NJ, Arntfield MD FRCPC R, Kory MD MPA P. *Point of Care Ultrasound.* 2nd ed. Elsevier; 2019:546.
6. Neskovic AN, Skinner H, Price S, et al. Focus cardiac ultrasound core curriculum and core syllabus of the European Association of Cardiovascular Imaging. *Eur Heart J Cardiovasc Imaging.* 2018;19(5):475-481. doi:10.1093/ehjci/jeu006
7. Núñez-Ramos JA, Pana-Tolozza MC, Palacio-Held SC. E-Point Septal Separation Accuracy for the Diagnosis of Mild and Severe Reduced Ejection Fraction in Emergency Department Patients. *POCUS.* 2022;7(1):160-165. doi:10.24908/pocus.v7i1.15220


Interactive Image Review



25




26




27



28




29




30




31




32



33




34




35




36



37




38




39




40




41




42



43




44



45



46



47

Please break into groups of 3 for the practical session

The slide features a black background with a white text instruction. At the top, there is a horizontal bar with a green segment on the left and a gray segment on the right. At the bottom, there is a green bar containing the ACP logo and tagline.

48